



DEGANIA™, a QMD™ product line trusted for its expertise in enteral feeding products, has introduced its latest advancement in Nasogastric Tubes. The newest design features high grade biocompatible radiopaque polyurethane material for easy insertion and patient comfort.

.....●

NASOGASTRIC FEEDING TUBES



PRODUCT DESCRIPTION

A Nasogastric Feeding Tube (NG) is a flexible device designated for patients who require intermittent or continuous feeding and/or medication delivery.

The DEGANIA™ line of Single-Port and Bi-Port Nasogastric Feeding Tubes are manufactured from premium materials designed to comply with ENFit® connection standards for improved patient safety.



Main Features:

- **Proximal Y-site connector** that allows for tube flushing and adding medication without disconnecting the feeding set
- **Large holes and open distal end** provide good fluid flow and minimizes clogging
- **Hydrophilic coating** for easy insertion and increase of patient discomfort
- **Radiopaque stripe** gives reference markings
- **Weighted/Non-Weighted, with and without stylet designs**
- **Pediatric & Adult range**

QMD™ offers private label program, saving lengthy and costly development, production and approval time. QMD™ Expertise labeled with your name.

Product Catalog

Reference Number	Product Description	Size (FR)	Length (mm)	Weight	Stylet	Sterile	Box Qty.	Carton Qty.
16605000580AQ	Nasogastric Feeding Tubes Single Port- Pediatric	5	500	No	No	Yes	25 Units	100 Units
16605030680AQ		6	500	No	Yes	Yes	25 Units	100 Units
16603030680AQ		6	900	No	Yes	Yes	25 Units	100 Units
16602030880AQ	Nasogastric Feeding Tubes Double Port	8	1,200	No	Yes	Yes	25 Units	100 Units
16601030880AQ		8	1,200	Yes	Yes	Yes	25 Units	100 Units
16602031080AQ		10	1,200	No	Yes	Yes	25 Units	100 Units
16601031080AQ		10	1,200	Yes	Yes	Yes	25 Units	100 Units
16602031280AQ		12	1,200	No	Yes	Yes	25 Units	100 Units
16601031280AQ		12	1,200	Yes	Yes	Yes	25 Units	100 Units
16602001280AQ		12	1,200	Yes	No	Yes	25 Units	100 Units
16602031480AQ		14	1,200	No	Yes	Yes	25 Units	100 Units
16602001480AQ		14	1,200	Yes	No	Yes	25 Units	100 Units

CONTACT US >>>



| www.qco.net/qmd



| qmeddevices@qco.net



| GLOBAL HEADQUARTERS
330100 Chagrin Blvd.
Suite 201
Pepper Pike, OH 44124
USA