



Degania™
Biometrix™
Arthesys™



SIGNAL CATHETER®

**Reduce premature balloon inflation risk
with a one-step secure insertion!**

Innovative catheterization technology
to help provide safe and efficient therapy.

The design reduces the risk of urethral
trauma from inadvertent inflation of the
catheter balloon in the urethra.



Degania™
DESIGNED TO DELIVER

The Signal Catheter® utilizes a proprietary, patented "Signal Balloon." This additional balloon provides several unique advantages over today's urinary catheters:



Serves as an "overflow" to limit the pressure in the primary catheter balloon if the catheter is improperly placed (ie. in the urethra).



When inflated, the "Signal Balloon" provides a visual indicator to the operator that the catheter is improperly placed and needs repositioning.



The catheter functions identically to current urinary catheters when properly placed, minimizing required operator training.

Signal Catheter® combines safety, performance, and ease of use.

MAIN FEATURES:

- **Signal Balloon** Limits urethral pressure to a safe level reducing the risk of urethral trauma.
- **100% silicone design** Suited for long-term urinary drainage.
- **Smooth, round shaft** Minimizes trauma during insertion and withdrawal.
- **Transparent shaft with radiopaque stripe** Enables urine drainage to be monitored and facilitates radiological imaging.
- **Soft shaft with closely fitting balloon**
Ensures easy catheter insertion and minimizes discomfort.
- **Thin, symmetrically inflated retention balloon**
Improves retention of the catheter in the bladder and prevents urine leakage.
- **Closed and curved radiopaque tip** Helps move through the curves of the urethra more easily.
- **Olive-shaped end** Allows for easy insertion around obstacles encountered in the urethra to minimize trauma.
- **Large cross-section of inner lumen** Enhances drainage performance.
- **Color-coded self-closing valve** Permits easy size indication.



QMD™ | Qure Medical Devices offers private label program, saving lengthy and costly development, production and approval time.

QMD Expertise labeled with your name.

CONTACT US >>>



Degania™
Biometrix™
Arthesys™

Degania™
DESIGNED TO DELIVER



| www.qmd.net



| marketing@qmd.net



GLOBAL HEADQUARTERS
30100 Chagrin Blvd.
Suite 201
Pepper Pike, OH 44124
USA

8390047 REV0