

Performing by design

Y Connector

The original vascular access system

Unique push pull system
Clean, safe and fast
Essential for Drug Eluting Stent




ARTHESYS

Committed to excellence in developing Arterial Therapeutic Systems

Y Connector

The original vascular access system



> A unique design

For easy, fast and accurate single handed intra vascular device insertion.

> Essential for DES

Perfect insertion device respecting the integrity of the coated surface.

> Clean and safe

Combination of haemostatic valve and integrated push-pull button valve opener. Cleaner and safer procedure with minimized backflow, limiting blood loss.



Technical data

Y-Connector

Inner Lumen Compatibility	7F / 0.092" / 2.33mm maximum
Pressure Resistance ⁽¹⁾ after introduction of a device	8 atm. maximum
Available range	20 & 50cm side arms No side arm Sterile or non sterile (in bulk)

⁽¹⁾ Haemostatic valve resistance measured in vitro, by applying injection pressure through the side port.

Accessories

Metallic Insertion tool ⁽²⁾

External Diameter	2.7F / 0.90mm
Inner Lumen Diameter	1.9F / 0.025" / 0.63mm
Usable length	10cm

Plastic Peelable Insertion Tool ⁽³⁾

Length	9cm
Inner Lumen Diameter	3.4F / 0.044" / 1.13mm

Torquer

Guidewire compatibility	From 0.010" to 0.035"
-------------------------	-----------------------

⁽²⁾ For the entire range

⁽³⁾ No side arm model only



20, rue Traversière - 92230 GENNEVILLIERS - FRANCE
Tel: + 33 (0)1 4111 8777 - info@arthesys.com

Take a look at our most up-to-date developments and latest news
<http://www.arthesys.com>